

For Research Use Only

BMP6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83269-2-PBS

Basic Information

Catalog Number: 83269-2-PBS	GenBank Accession Number: NM_001718	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 654	CloneNo.: 240109F7
Source: Rabbit	UNIPROT ID: P22004	Affinity: $K_D = 1.54 \times 10^{-9}$
Isotype: IgG	Full Name: bone morphogenetic protein 6	$K_{off} = 1.28 \times 10^{-3}$
	Calculated MW: 57 kDa	$K_{on} = 8.28 \times 10^5$

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C .
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00226-1, BMP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83269-2-PBS. Detection antibody: 83269-3-PBS. Standard: SY01500. Range: 0.625-20 ng/mL