

For Research Use Only

# Cytochrome c Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83276-4-PBS



## Basic Information

Catalog Number:

83276-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1455

GenBank Accession Number:

BC009578

GeneID (NCBI):

54205

UNIPROT ID:

P99999

Full Name:

cytochrome c, somatic

Calculated MW:

12 kDa

Purification Method:

Protein A purification

CloneNo.:

230503G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01934-1, Cytochrome c Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83276-4-PBS. Detection antibody: 83276-5-PBS. Standard: Ag1455. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01934-2, Cytochrome c Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83276-3-PBS. Detection antibody: 83276-4-PBS. Standard: Ag1455. Range: 0.781-100 ng/mL.