

For Research Use Only

# CCL17 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83278-1-PBS

## Basic Information

Catalog Number: 83278-1-PBS	GenBank Accession Number: BC112066	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 6361	CloneNo.: 230354E8
Source: Rabbit	UNIPROT ID: Q92583	Affinity: $K_D = 7.69 \times 10^{-11}$
Isotype: IgG	Full Name: chemokine (C-C motif) ligand 17	$K_{off} = 1.38 \times 10^{-4}$ $K_{on} = 1.80 \times 10^6$
	Calculated MW: 94 aa, 11 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00308-1, CCL17 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83278-1-PBS. Detection antibody: 83278-2-PBS. Standard: Eg0379. Range: 0.625-20 ng/mL