

For Research Use Only

CCL17/TARC Recombinant antibody, PBS Only (Detector)

Catalog Number: 83278-2-PBS



Basic Information

Catalog Number:

83278-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC112066

GeneID (NCBI):

6361

UNIPROT ID:

Q92583

Full Name:

chemokine (C-C motif) ligand 17

Calculated MW:

94 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

230354C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00308-1, CCL17 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83278-1-PBS. Detection antibody: 83278-2-PBS. Standard: Eg0379. Range: 0.625-20 ng/mL