

For Research Use Only

Mouse Albumin Recombinant antibody, PBS Only (Detector)



Catalog Number: 83289-2-PBS

Basic Information

Catalog Number: 83289-2-PBS	GenBank Accession Number: NM_009654.4	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 11657	CloneNo.: 240245H4
Source: Rabbit	UNIPROT ID: P07724	
Isotype: IgG	Full Name: albumin	
	Calculated MW: 69 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00329-2, MOUSE Albumin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83289-1-PBS. Detection antibody: 83289-2-PBS. Standard: Eg0886. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00329-3, MOUSE Albumin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83289-3-PBS. Detection antibody: 83289-2-PBS. Standard: Eg0886. Range: 0.313-40 ng/mL.