

For Research Use Only

Mouse TREM-1/CD354 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83294-1-PBS



Basic Information

Catalog Number:

83294-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0874

GenBank Accession Number:

NM_021406.5

GeneID (NCBI):

58217

UNIPROT ID:

Q9JKE2

Full Name:

triggering receptor expressed on myeloid cells 1

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

240219D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00316-1, MOUSE Trem1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83294-1-PBS. Detection antibody: 83294-2-PBS. Standard: Eg0874. Range: 1.25-40 ng/mL