

For Research Use Only

Mouse CD14 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83297-3-PBS

Basic Information

Catalog Number:

83297-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009841.4

GeneID (NCBI):

12475

UNIPROT ID:

P10810

Full Name:

CD14 antigen

Calculated MW:

39 kDa

Purification Method:

Protein A purification

CloneNo.:

240225E7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00313-1, MOUSE CD14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83297-3-PBS. Detection antibody: 83297-1-PBS. Standard: Eg0847. Range: 1.25-40 ng/mL