

For Research Use Only

STX1B Recombinant antibody, PBS Only (Capture)



Catalog Number: 83298-5-PBS

Basic Information

Catalog Number: 83298-5-PBS	GenBank Accession Number: BC062298	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 112755	CloneNo.: 240156G5
Source: Rabbit	UNIPROT ID: P61266	
Isotype: IgG	Full Name: syntaxin 1B	
Immunogen Catalog Number: AG10275	Calculated MW: 288 aa, 33 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00312-2, STX1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83298-5-PBS. Detection antibody: 83298-4-PBS. Standard: Ag10275. Range: 1.25-40 ng/mL.