

For Research Use Only

TNFSF18 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83299-2-PBS



Basic Information

Catalog Number:

83299-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20979

GenBank Accession Number:

BC069319

GeneID (NCBI):

8995

UNIPROT ID:

Q9UNG2

Full Name:

tumor necrosis factor (ligand)

superfamily, member 18

Calculated MW:

177 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

240098G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01084-4, TNFSF18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83299-1-PBS. Detection antibody: 83299-2-PBS. Standard: Ag20979. Range: 3.125-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83299-2-PBS against Human TNFSF18 was performed. The affinity constant is 1.02 nM.