

For Research Use Only

SMARCA4/BRG1 Recombinant antibody, PBS Only

Catalog Number: 83310-4-PBS



Basic Information

Catalog Number:

83310-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16256

GenBank Accession Number:

BC150298

GeneID (NCBI):

6597

UNIPROT ID:

P51532

Full Name:

SWI/SNF related, matrix associated,
actin dependent regulator of
chromatin, subfamily a, member 4

Calculated MW:

1647 aa, 185 kDa

Purification Method:

Protein A purification

CloneNo.:

240172H1

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

