

For Research Use Only

SMARCA4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83310-5-PBS

Basic Information

Catalog Number:

83310-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16256

GenBank Accession Number:

BC150298

GeneID (NCBI):

6597

UNIPROT ID:

P51532

Full Name:

SWI/SNF related, matrix associated,
actin dependent regulator of
chromatin, subfamily a, member 4

Calculated MW:

1647 aa, 185 kDa

Purification Method:

Protein A purification

CloneNo.:

240172E7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00334-2, SMARCA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83310-1-PBS. Detection antibody: 83310-5-PBS. Standard: Ag16256. Range: 0.625-20 ng/mL.