

For Research Use Only

# BRAP Recombinant antibody, PBS Only (Capture)

Catalog Number: 83314-3-PBS



## Basic Information

Catalog Number:

83314-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26728

GenBank Accession Number:

BC136698

GeneID (NCBI):

8315

UNIPROT ID:

Q7Z569

Full Name:

BRCA1 associated protein

Calculated MW:

592 aa, 67 kDa

Purification Method:

Protein A purification

CloneNo.:

240099A7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00338-2, BRAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83314-3-PBS. Detection antibody: 83314-4-PBS. Standard: Ag26728. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00338-3, BRAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83314-3-PBS. Detection antibody: 83314-1-PBS. Standard: Ag26728. Range: 0.313-40 ng/mL