

For Research Use Only

TRIM34 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83316-4-PBS

Basic Information

| | | |
|--------------------------------------|--|--|
| Catalog Number: 83316-4-PBS | GenBank Accession Number: NM_001003827 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 53840 | CloneNo.: 240180F12 |
| Source: Rabbit | UNIPROT ID: Q9BYJ4 | |
| Isotype: IgG | Full Name: tripartite motif-containing 34 | |
| Immunogen Catalog Number: AG34444 | Calculated MW: 57 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00340-2, TRIM34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83316-4-PBS. Detection antibody: 83316-3-PBS. Standard: Ag34444. Range: 1.25-40 ng/mL.