

For Research Use Only

# ADRB2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83317-4-PBS

## Basic Information

Catalog Number: 83317-4-PBS	GenBank Accession Number: BC012481	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 154	CloneNo.: 240186D9
Source: Rabbit	UNIPROT ID: P07550	
Isotype: IgG	Full Name: adrenergic, beta-2-, receptor, surface	
Immunogen Catalog Number: AG3715	Calculated MW: 46 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00341-3, ADRB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83317-4-PBS. Detection antibody: 83317-3-PBS. Standard: Ag3715. Range: 0.156-20 ng/mL