

For Research Use Only

APIP Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 83320-3-PBS



Basic Information

Catalog Number:

83320-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2729

GenBank Accession Number:

BC009077

GeneID (NCBI):

51074

UNIPROT ID:

Q96GX9

Full Name:

APAF1 interacting protein

Calculated MW:

242 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

240189F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00346-1, APDP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83320-3-PBS. Detection antibody: 83320-1-PBS. Standard: Ag2729. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00346-3, APDP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83320-3-PBS. Detection antibody: 83320-2-PBS. Standard: Ag2729. Range: 0.313-40 ng/mL