

For Research Use Only

APIP Recombinant antibody

Catalog Number: 83320-5-RR



Basic Information

Catalog Number:

83320-5-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2729

GenBank Accession Number:

BC009077

GeneID (NCBI):

51074

UNIPROT ID:

Q96GX9

Full Name:

APAF1 interacting protein

Calculated MW:

242 aa, 27 kDa

Observed MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

240189E2

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : PC-3 cells, K-562 cells.

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83320-5-RR (AP-1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 83320-5-RR against Human AP-1 were performed. The affinity constant is 5.86 nM.