

For Research Use Only

# EGLN3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83321-4-PBS



## Basic Information

Catalog Number:

83321-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27464

GenBank Accession Number:

BC010992

GeneID (NCBI):

112399

UNIPROT ID:

Q9H6Z9

Full Name:

egl nine homolog 3 (C. elegans)

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

240165A12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00347-1, EGLN3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83321-1-PBS. Detection antibody: 83321-4-PBS. Standard: Ag27464. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83321-4-PBS against Human EGLN3 was performed. The affinity constant is below 1 pM.