

For Research Use Only

Bcl-XL Recombinant antibody, PBS Only (Detector)

Catalog Number: 83330-5-PBS



Basic Information

Catalog Number:

83330-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25604

GenBank Accession Number:

BC019307

GeneID (NCBI):

598

UNIPROT ID:

Q07817

Full Name:

BCL2-like 1

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

240450H1

Applications

Tested Applications:

WB, Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00460-2, Bcl-xL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83330-4-PBS. Detection antibody: 83330-5-PBS. Standard: Ag25604. Range: 0.625-80 ng/mL