

For Research Use Only

CMTM6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83334-1-PBS

Basic Information

Catalog Number: 83334-1-PBS	GenBank Accession Number: BC002797	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54918	CloneNo.: 240194A12
Source: Rabbit	UNIPROT ID: Q9NX76	
Isotype: IgG	Full Name: CKLF-like MARVEL transmembrane domain containing 6	
Immunogen Catalog Number: AG34445		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00352-1, CMTM6 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83334-2-PBS. Detection antibody: 83334-1-PBS. Standard: Ag34445. Range: 2.5-80 ng/mL.