

For Research Use Only

Mouse IL-1RA Recombinant antibody, PBS Only (Detector)

Catalog Number: 83335-1-PBS



Basic Information

Catalog Number:

83335-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_031167.5

GeneID (NCBI):

16181

UNIPROT ID:

Q542W1

Full Name:

interleukin 1 receptor antagonist

Calculated MW:

18KD

Purification Method:

Protein A purification

CloneNo.:

240208A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

The interleukin-1 receptor antagonist (IL-1RN) belongs to the IL-1 family. L1RN is a protein that binds to IL1 receptors, and inhibits the binding of IL1-alpha and IL1-beta.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00353-1, MOUSE IL-1RA Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83335-2-PBS. Detection antibody: 83335-1-PBS. Standard: Eg0685. Range: 0.625-80 ng/mL.