

For Research Use Only

DNASE1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83337-4-PBS

Basic Information

Catalog Number: 83337-4-PBS	GenBank Accession Number: BC029437	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1773	CloneNo.: 240222A8
Source: Rabbit	UNIPROT ID: P24855	
Isotype: IgG	Full Name: deoxyribonuclease I	
Immunogen Catalog Number: AG10117	Calculated MW: 282 aa, 31 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00356-2, DNASE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83337-2-PBS. Detection antibody: 83337-4-PBS. Standard: Ag10117. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00356-3, DNASE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83337-3-PBS. Detection antibody: 83337-4-PBS. Standard: Ag10117. Range: 0.625-40 ng/mL.