

For Research Use Only

DNASE1 Recombinant antibody, PBS Only

Catalog Number: 83337-6-PBS



Basic Information

Catalog Number:

83337-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10117

GenBank Accession Number:

BC029437

GeneID (NCBI):

1773

UNIPROT ID:

P24855

Full Name:

deoxyribonuclease I

Calculated MW:

282 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

240222F11

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 83337-6-RR (DNase I antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83337-6-PBS in a different storage buffer formulation.