

For Research Use Only

# CDH9 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83338-2-PBS

## Basic Information

Catalog Number: 83338-2-PBS	GenBank Accession Number: BC113745	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1007	CloneNo.: 240289F10
Source: Rabbit	UNIPROT ID: Q9ULB4	
Isotype: IgG	Full Name: cadherin 9, type 2 (T1-cadherin)	
Immunogen Catalog Number: AG16437	Calculated MW: 789 aa, 89 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00359-1, CDH9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83338-1-PBS. Detection antibody: 83338-2-PBS. Standard: Ag16437. Range: 0.313-40 ng/mL