

For Research Use Only

ST3GAL5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83339-1-PBS

Basic Information

Catalog Number:

83339-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6207

GenBank Accession Number:

BC065936

GeneID (NCBI):

8869

UNIPROT ID:

Q9UNP4

Full Name:

ST3 beta-galactoside alpha-2,3-sialyltransferase 5

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

240148A9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00381-2, ST3GAL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83339-1-PBS. Detection antibody: 83339-3-PBS. Standard: Ag6207. Range: 0.625-40 ng/mL