

For Research Use Only

# ST3GAL5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83339-3-PBS

## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
83339-3-PBS	BC065936	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	8869	240148C10
Source:	UNIPROT ID:	
Rabbit	Q9UNP4	
Isotype:	Full Name:	
IgG	ST3 beta-galactoside alpha-2,3-sialyltransferase 5	
Immunogen Catalog Number:	Calculated MW:	
AG6207	48 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00381-1, ST3GAL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83339-2-PBS. Detection antibody: 83339-3-PBS. Standard: Ag6207. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00381-2, ST3GAL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83339-1-PBS. Detection antibody: 83339-3-PBS. Standard: Ag6207. Range: 0.625-40 ng/mL