

For Research Use Only

CRIP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83340-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
83340-1-PBS	BC002738	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	1396	240176C5
Source:	UNIPROT ID:	
Rabbit	P50238	
Isotype:	Full Name:	
IgG	cysteine-rich protein 1 (intestinal)	
Immunogen Catalog Number:	Calculated MW:	
AG7588	9 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00382-2, CRIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83340-1-PBS. Detection antibody: 83340-3-PBS. Standard: Ag7588. Range: 0.156-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83340-1-PBS against Human CRIP1 was performed. The affinity constant is 1.32 nM.