

For Research Use Only

C1QTNF2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83343-2-PBS

Basic Information

Catalog Number:

83343-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9183

GenBank Accession Number:

BC011699

GeneID (NCBI):

114898

UNIPROT ID:

Q9BXJ5

Full Name:

C1q and tumor necrosis factor related
protein 2

Calculated MW:

330 aa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

240330F8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00385-1, C1QTNF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83343-3-PBS. Detection antibody: 83343-2-PBS. Standard: Ag9183. Range: 0.313-40 ng/mL