

For Research Use Only

# C1QTNF2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83343-4-PBS



## Basic Information

Catalog Number:

83343-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9183

GenBank Accession Number:

BC011699

GeneID (NCBI):

114898

UNIPROT ID:

Q9BXJ5

Full Name:

C1q and tumor necrosis factor related protein 2

Calculated MW:

330 aa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

240330G2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00385-2, C1QTNF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83343-1-PBS. Detection antibody: 83343-4-PBS. Standard: Ag9183. Range: 0.625-40 ng/mL