

For Research Use Only

MTF1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83344-2-PBS

Basic Information

| | | | | | |
|---------------------------|-------------|---------------------------|---------------------------------------|----------------------|------------------------|
| Catalog Number: | 83344-2-PBS | GenBank Accession Number: | BC014454 | Purification Method: | Protein A purification |
| Size: | 1 mg/ml | GeneID (NCBI): | 4520 | CloneNo.: | 240336H1 |
| Source: | Rabbit | UNIPROT ID: | Q14872 | | |
| Isotype: | IgG | Full Name: | metal-regulatory transcription factor | | |
| Immunogen Catalog Number: | AG21848 | | 1 | | |
| | | Calculated MW: | 753 aa, 81 kDa | | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00386-2, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-2-PBS. Detection antibody: 83344-3-PBS. Standard: Ag21848. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00386-3, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-2-PBS. Detection antibody: 83344-1-PBS. Standard: Ag21848. Range: 0.625-20 ng/mL