

For Research Use Only

MTF1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83344-3-PBS

Basic Information

Catalog Number:	83344-3-PBS	GenBank Accession Number:	BC014454	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	4520	CloneNo.:	240336F7
Source:	Rabbit	UNIPROT ID:	Q14872		
Isotype:	IgG	Full Name:	metal-regulatory transcription factor		
Immunogen Catalog Number:	AG21848		1		
		Calculated MW:	753 aa, 81 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00386-1, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-1-PBS. Detection antibody: 83344-3-PBS. Standard: Ag21848. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00386-2, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-2-PBS. Detection antibody: 83344-3-PBS. Standard: Ag21848. Range: 0.313-20 ng/mL