

For Research Use Only

LYPD3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83345-2-PBS

Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 83345-2-PBS | GenBank Accession Number: NM_014400 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 27076 | CloneNo.: 240347F4 |
| Source: Rabbit | UNIPROT ID: O95274 | |
| Isotype: IgG | Full Name: LY6/PLAUR domain containing 3 | |
| Immunogen Catalog Number: AG34961 | Calculated MW: 36 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00387-1, LYPD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83345-3-PBS. Detection antibody: 83345-2-PBS. Standard: Ag34961. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00387-2, LYPD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83345-1-PBS. Detection antibody: 83345-2-PBS. Standard: Ag34961. Range: 0.156-20 ng/mL