

For Research Use Only

# LYPD3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83345-3-PBS

## Basic Information

Catalog Number: 83345-3-PBS	GenBank Accession Number: NM_014400	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 27076	CloneNo.: 240347D8
Source: Rabbit	UNIPROT ID: O95274	
Isotype: IgG	Full Name: LY6/PLAUR domain containing 3	
Immunogen Catalog Number: AG34961	Calculated MW: 36 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00387-1, LYPD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83345-3-PBS. Detection antibody: 83345-2-PBS. Standard: Ag34961. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00387-3, LYPD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83345-1-PBS. Detection antibody: 83345-3-PBS. Standard: Ag34961. Range: 0.156-20 ng/mL