

For Research Use Only

SLC38A7 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83346-2-PBS

Basic Information

Catalog Number: 83346-2-PBS	GenBank Accession Number: BC001961	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55238	CloneNo.: 240147G5
Source: Rabbit	UNIPROT ID: Q9NVC3	
Isotype: IgG	Full Name: solute carrier family 38, member 7	
Immunogen Catalog Number: AG35081	Calculated MW: 50 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00351-1, SLC38A7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83346-1-PBS. Detection antibody: 83346-2-PBS. Standard: Ag35081. Range: 0.313-40 ng/mL.