

For Research Use Only

# TCIRG1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83351-1-PBS

## Basic Information

Catalog Number:

83351-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3378

GenBank Accession Number:

BC018133

GeneID (NCBI):

10312

UNIPROT ID:

Q13488

Full Name:

T-cell, immune regulator 1, ATPase,  
H<sup>+</sup> transporting, lysosomal VO subunit  
A3

Calculated MW:

830 aa, 93 kDa

Purification Method:

Protein A purification

CloneNo.:

240166A4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00365-3, TCIRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83351-1-PBS. Detection antibody: 83351-2-PBS. Standard: Ag3378. Range: 0.313-40 ng/mL