

For Research Use Only

TCIRG1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83351-5-PBS

Basic Information

Catalog Number:

83351-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3378

GenBank Accession Number:

BC018133

GeneID (NCBI):

10312

UNIPROT ID:

Q13488

Full Name:

T-cell, immune regulator 1, ATPase,
H⁺ transporting, lysosomal VO subunit
A3

Calculated MW:

830 aa, 93 kDa

Purification Method:

Protein A purification

CloneNo.:

240166F7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00365-1, TCIRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83351-3-PBS. Detection antibody: 83351-5-PBS. Standard: Ag3378. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00365-2, TCIRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83351-4-PBS. Detection antibody: 83351-5-PBS. Standard: Ag3378. Range: 0.313-40 ng/mL