

For Research Use Only

# CNKS2 Recombinant antibody, PBS Only (Capture)



Catalog Number:83354-2-PBS

## Basic Information

Catalog Number:

83354-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21174

GenBank Accession Number:

BC126121

GeneID (NCBI):

22866

UNIPROT ID:

Q8WXL2

Full Name:

connector enhancer of kinase

suppressor of Ras 2

Calculated MW:

1034 aa, 118 kDa

Purification Method:

Protein A purification

CloneNo.:

240350E2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00389-1, CNKSR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83354-2-PBS. Detection antibody: 83354-3-PBS. Standard: Ag21174. Range: 3.125-100ng/mL