

For Research Use Only

CNKS2 Recombinant antibody, PBS Only (Detector)



Catalog Number:83354-4-PBS

Basic Information

Catalog Number:	83354-4-PBS	GenBank Accession Number:	BC126121	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	22866	CloneNo.:	240350H8
Source:	Rabbit	UNIPROT ID:	Q8WXL2		
Isotype:	IgG	Full Name:	connector enhancer of kinase suppressor of Ras 2		
Immunogen Catalog Number:	AG21174	Calculated MW:	1034 aa, 118 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00389-2, CNKSR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83354-1-PBS. Detection antibody: 83354-4-PBS. Standard: Ag21174. Range: 6.25-100ng/mL.