

For Research Use Only

DUSP6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83355-2-PBS



Basic Information

Catalog Number: 83355-2-PBS	GenBank Accession Number: BC003562	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1848	CloneNo.: 240303C4
Source: Rabbit	UNIPROT ID: Q16828	
Isotype: IgG	Full Name: dual specificity phosphatase 6	
Immunogen Catalog Number: AG0705	Calculated MW: 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00391-1, DUSP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83355-1-PBS. Detection antibody: 83355-2-PBS. Standard: Ag0705. Range: 2.5-80ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83355-2-PBS against Human DUSP6 was performed. The affinity constant is 1.87 nM.