

For Research Use Only

DUSP6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83356-4-PBS



Basic Information

Catalog Number:

83356-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23474

GenBank Accession Number:

BC003562

GeneID (NCBI):

1848

UNIPROT ID:

Q16828

Full Name:

dual specificity phosphatase 6

Calculated MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

240295G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00391-4, DUSP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83356-3-PBS. Detection antibody: 83356-4-PBS. Standard: Ag23474. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00391-6, DUSP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83356-2-PBS. Detection antibody: 83356-4-PBS. Standard: Ag23474. Range: 0.625-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83356-4-PBS against Human DUSP6 was performed. The affinity constant is 0.834 nM.