

For Research Use Only

SLC25A46 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83360-3-PBS

Basic Information

Catalog Number:

83360-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27293

GenBank Accession Number:

BC017169

GeneID (NCBI):

91137

UNIPROT ID:

Q96AG3

Full Name:

solute carrier family 25, member 46

Calculated MW:

418 aa, 46 kDa

Purification Method:

Protein A purification

CloneNo.:

240114A9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00369-1, SLC25A46 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83360-2-PBS. Detection antibody: 83360-3-PBS. Standard: Ag27293. Range: 0.625-80ng/mL.



Cytometric bead array standard curve of MP00369-2, SLC25A46 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83360-1-PBS. Detection antibody: 83360-3-PBS. Standard: Ag27293. Range: 2.5-80ng/mL.