

For Research Use Only

# NEDD9 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83362-1-PBS

## Basic Information

Catalog Number: 83362-1-PBS	GenBank Accession Number: BC040207	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4739	CloneNo.: 240310A2
Source: Rabbit	UNIPROT ID: Q14511	
Isotype: IgG	Full Name: neural precursor cell expressed, developmentally down-regulated 9	
Immunogen Catalog Number: AG34170	Calculated MW: 834 aa, 93 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00378-1, NEDD9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83362-3-PBS. Detection antibody: 83362-1-PBS. Standard: Ag34170. Range: 0.625-80 ng/mL