

For Research Use Only

# PCARE Recombinant antibody, PBS Only (Capture)



Catalog Number: 83364-3-PBS

## Basic Information

Catalog Number: 83364-3-PBS	GenBank Accession Number: NM_001029883	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 388939	CloneNo.: 240352B6
Source: Rabbit	UNIPROT ID: A6NGG8	
Isotype: IgG	Full Name: chromosome 2 open reading frame 71	
Immunogen Catalog Number: AG34774	Calculated MW: 140 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00392-1, PCARE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83364-3-PBS. Detection antibody: 83364-1-PBS. Standard: Ag34774. Range: 1.25-40 ng/mL.