

For Research Use Only

# Mouse VEGFR2/KDR Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83366-2-PBS



## Basic Information

Catalog Number:

83366-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_010612.2

GeneID (NCBI):

16542

UNIPROT ID:

P35918-1

Full Name:

kinase insert domain protein receptor

Calculated MW:

150KD

Purification Method:

Protein A purification

CloneNo.:

240214G8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00394-2, MOUSE VEGFR2/KDR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83366-2-PBS. Detection antibody: 83366-3-PBS. Standard: Eg0703. Range: 2.5-80 ng/mL.



Cytometric bead array standard curve of MP00394-3, MOUSE VEGFR2/KDR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83366-1-PBS. Detection antibody: 83366-2-PBS. Standard: Eg0703. Range: 5-80 ng/mL.