

For Research Use Only

Mouse Afp Recombinant antibody, PBS Only (Detector)



Catalog Number: 83367-3-PBS

Basic Information

Catalog Number:

83367-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_007423.4

GeneID (NCBI):

11576

UNIPROT ID:

P02772

Full Name:

alpha fetoprotein

Calculated MW:

67KD

Purification Method:

Protein A purification

CloneNo.:

240232B1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00395-3, MOUSE AFP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83367-5-PBS. Detection antibody: 83367-3-PBS. Standard: Eg0670. Range: 2.5-40 ng/mL.