

For Research Use Only

Rat MCP-1/CCL2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83370-3-PBS



Basic Information

Catalog Number:

83370-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_031530

GeneID (NCBI):

24770

UNIPROT ID:

P14844

Full Name:

chemokine (C-C motif) ligand 2

Calculated MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

240341E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00390-1, RAT MCP-1/CCL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83370-3-PBS. Detection antibody: 83370-1-PBS. Standard: Eg0397. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00390-2, RAT MCP-1/CCL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83370-3-PBS. Detection antibody: 83370-2-PBS. Standard: Eg0397. Range: 0.156-20 ng/mL