

For Research Use Only

BAP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83371-1-PBS



Basic Information

Catalog Number:

83371-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34415

GenBank Accession Number:

BC001596

GeneID (NCBI):

8314

UNIPROT ID:

Q92560

Full Name:

BRCA1 associated protein-1 (ubiquitin carboxy-terminal hydrolase)

Calculated MW:

81-91 kDa

Purification Method:

Protein A purification

CloneNo.:

240349A7

Applications

Tested Applications:

WB, Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.



Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83371-1-RR (BAP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83371-1-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP00357-1, BAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83371-1-PBS. Detection antibody: 83371-2-PBS. Standard: Ag34415. Range: 1.25-80 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 83371-1-RR against Human BAP1 were performed. The affinity constant is below 1 pM.