

For Research Use Only

BAP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83371-2-PBS



Basic Information

Catalog Number:

83371-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34415

GenBank Accession Number:

BC001596

GeneID (NCBI):

8314

UNIPROT ID:

Q92560

Full Name:

BRCA1 associated protein-1 (ubiquitin carboxy-terminal hydrolase)

Calculated MW:

81-91 kDa

Purification Method:

Protein A purification

CloneNo.:

240349F5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00357-1, BAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83371-1-PBS. Detection antibody: 83371-2-PBS. Standard: Ag34415. Range: 1.25-80 ng/mL.