

For Research Use Only

KIAA1199 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83373-2-PBS



Basic Information

Catalog Number:

83373-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15527

GenBank Accession Number:

BC020256

GeneID (NCBI):

57214

UNIPROT ID:

Q8WUJ3

Full Name:

KIAA1199 / CEMIP

Calculated MW:

1361 aa, 153 kDa

Purification Method:

Protein A purification

CloneNo.:

240372F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00373-2, KIAA1199 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83373-1-PBS. Detection antibody: 83373-2-PBS. Standard: Ag15527. Range: 1.25-80 ng/mL