

For Research Use Only

Mouse ICOSLG/B7-H2/CD275 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83375-3-PBS



Basic Information

Catalog Number:

83375-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_015790.3

GeneID (NCBI):

50723

UNIPROT ID:

Q9JHJ8-1

Full Name:

icos ligand

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

240243D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00376-1, MOUSE ICOSLG/B7-H2/CD275 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83375-2-PBS. Detection antibody: 83375-3-PBS. Standard: Eg0923. Range: 0.313-40 ng/mL