

For Research Use Only

# Mouse ICOSLG/B7-H2/CD275 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83375-4-PBS



## Basic Information

Catalog Number:

83375-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_015790.3

GeneID (NCBI):

50723

UNIPROT ID:

Q9JHJ8-1

Full Name:

icos ligand

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

240243H6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00376-3, MOUSE ICOSLG/B7-H2/CD275 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83375-1-PBS. Detection antibody: 83375-4-PBS. Standard: Eg0923. Range: 0.625-40 ng/mL